


Concentrated Washroom Cleaner

Eco Select™ Concentrated Washroom Cleaner provides control to routine hygiene operations. The specially designed bottle enables exact dosing and optimised economy. Kills 99.99% of germs including Listeria, MRSA and E.Coli*.

The fresh fragranced formulation cuts through grease, grime and soap scum and removes limescale from sinks and taps leaving them gleaming and fresh. When diluted 1:60 it provides cleaning and disinfection of all hard surfaces including tiles, sinks, toilet bowls and urinals.

- Kills 99.99% of germs (passes BS EN 1276)
- Removes limescale & soap scum
- Cleans, disinfects & leaves a sparkling shine
- Forest fragrance
- 1Ltr makes 100 Trigger Sprays or 50 Buckets

(* Diluted 1:60).



1 Litre Makes

100
Trigger
Sprays



OR



50
Buckets

QUALITY ASSURANCE:

This product is manufactured in the UK by MIRIUS, A Coventry Group Company.

Produced under ISO 9001 Quality Management System & ISO 14001 Environmental Management System. This ensures our products and services are of the highest possible standard.

This product has not been tested on animals.

Date: 15/10/2019

BIODEGRADABILITY:

All surfactants used in Super Professional products comply with the current European regulations concerning biodegradability & protection of the environment.

- BOTTLE - HDPE**
widely recycled
- LABEL - PP**
widely recycled



Scan Me



Eco Select™ Concentrated Washroom provides control to routine hygiene operations, the specially designed bottle enables exact dosing and optimised economy, kills 99.99% of germs including Listeria, MRSA and E.Coli*.

The fresh fragranced formulation cuts through grease, grime and soap scum and removes limescale from sinks and taps leaving them gleaming and fresh. When diluted 1:60 it provides cleaning and disinfection of all hard surfaces including tiles, sinks, toilet bowls and urinals.

Cap Instructions: To remove cap: Place bottle on a firm flat surface. Push down cap and turn anti-clockwise. To replace cap: Place bottle on a firm flat surface. Screw cap clockwise until tight. If correctly closed, the cap will give a click if turned anti-clockwise.

Directions for use: Once open, squeeze bottle until measure fills to desired dose.

For trigger bottles: Add 600ml of clean water into the trigger bottle. Then add 1 x 10ml dose of product. Shake gently to mix. To ensure disinfection, leave for 5 minutes.

For general cleaning: Add 5 litres of warm water into a bucket. Then add 1 x 20ml dose of product. Stir to mix. Note: Use only as instructed. Do not use undiluted. * Diluted 1:60

Storage: Store sealed in original container in a cool, dry place away from sunlight.

Product not intended for animal use. Safety Data Sheet available on request.

For Safety Data Sheet please scan me:

BOTTLE - HDPE widely recycled
LABEL - PP widely recycled

Contains: Citric Acid Anhydrous, 2-Hal (C12-18) Dimethyl Benzyl Ammonium Chloride (ADBAC/BAC (C12-18)) 3.35g/100g, 2-(4-Hydroxyphenyl)Ethanol, Alcohols, C12-14, Ethylene Glycol.

This product has not been tested on animals.

1ltr e

Order Code: 800-304-0020

mirius
Woodhams Road, Siskin Drive, Coventry, CV3 4FX, Customer Care Line: +44(0)2476 639739 www.mirius.com


200-304-0020

ORDER CODE:

1ltr - 800-304-0020



Box qty per case:
4 x 1Ltr Bottles



140 cases to a pallet
5 layers to a pallet